


veterinary technology for life

Eickemeyer®

www.eickemeyer.com

Telephone +1 519 273 5558

Consumables



EICKWRAP BANDAGES

These high quality bandages have a self-textured surface structure allowing the layers of the wrap to bond together for maximum holding power. They do not constrict or impede the blood circulation and facilitate natural movement, whilst still providing support. EickWrap bandages are ideal for horses legs, hooves or dog and cat paws.

- ▶ Cohesive bandage
- ▶ Conforms itself to the shape of the area where it is applied
- ▶ Self-adhesive, but it will not stick to the skin
- ▶ Allows air to circulate over the skin and around the dressing
- ▶ Helps keep dressing in place
- ▶ Individually packaged – assorted colours

5 cm / 36 per box

4808031

7.5 cm / 24 per box

4808032

10 cm / 18 per box


4808033



4808031

4808032

4808033



GENIA STRAP BAND

- ▶ Cohesive elastic bandages with bitter taste
- ▶ Self-adhesive, without sticking to the skin or the hair
- ▶ Easy to cut, no scissors needed

▶ Dimensions: W 5 mm x L 4.5 m

▶ Pack of 24

4809476

▶ Dimensions: W 7.5 mm x L 4.5 m

▶ Pack of 12

4809477

▶ Dimensions: W 10 mm x L 4.5 m

▶ Pack of 12

4809471



4809476

GUTTER SPLINTS

- ▶ For temporary immobilisation of fractures
- ▶ Strong and lightweight plastic, green
- ▶ Comfortable, non-irritating foam padding
- ▶ Anatomically shaped, fully radiolucent
- ▶ Scoop end to be placed at toes
- ▶ Bag of 12

Dimensions (in cm): L 20 x W 3

4808070

Dimensions (in cm): L 30 x W 4

4808071

Dimensions (in cm): L 40 x W 5

4808072



ROLL SPLINTS

Aluminium roll splints with foam on both sides

- ▶ Easy to cut with serrated Mayo scissors to fit the animal's limb
- ▶ Radiolucent
- ▶ Flexible, atraumatic, easily takes the shape of the limb requiring immobilization
- ▶ Available in 3 widths
- ▶ Supplied in pack of four

Dimensions (in cm): W 4.5 x L 25

4804961

Dimensions (in cm): W 5.5 x L 35

4804962

Dimensions (in cm): W 9 x L 40

4804963



ALUMINIUM SPLINTS

These lightweight aluminium splints are laminated with sturdy foam and can be cut/formed to correctly immobilise most long bone injuries.

- ▶ Dimensions (in cm): W 1.5 x L 50
- ▶ Pack of 12
- ▶ Dimensions (in cm): W 2.5 x L 50
- ▶ Pack of 4

4804975

4804977

- ▶ Dimensions (in cm): W 2 x L 50
- ▶ Pack of 4

4804976



4804975

ABDOMINAL BANDAGE

- ▶ Soft and non abrasive on the skin
- ▶ 100 % cotton
- ▶ Washable at 95 °C
- ▶ Adjustable Velcro strap to ensure a perfect fit
- ▶ Tailored fit with 7 different sizes

Size 1

- ▶ Suitable for cats
- ▶ Length: 32 cm

4802764

Size 2

- ▶ Suitable for Poodles
- ▶ Length: 38 cm

4802765

Size 3

- ▶ Suitable for Setters
- ▶ Length: 47 cm

4802766

Size 4

- ▶ Suitable for Labradors
- ▶ Length: 51 cm

4802767

Size 5

- ▶ Suitable for German Shepherds
- ▶ Length: 60 cm

4802768

Size 6

- ▶ Suitable for Great Danes
- ▶ Length: 67 cm
- ▶ Width: 25 cm

4802769

Size 7

- ▶ Suitable for St. Bernards
- ▶ Length: 67 cm
- ▶ Width: 25 cm

4802770



4802769

SUGI SWAB

- ▶ Sugi absorbent materials are especially suitable for dabbing those areas where special care must be taken (e. g. neural structures and ophthalmic tissue)
- ▶ Sugi absorbent materials are manufactured from cellulose material
- ▶ The material undergoes soft/elastic expansion immediately upon contacting moisture
- ▶ The absorptive capacity is approximately 17 times its own weight
- ▶ Wedge tip
- ▶ 7 mm
- ▶ Set of 200 (Packs of 10 – sterile)

4801417



VACUUM IMOBILIZATION MATTRESS


- ▶ Radiolucent, easy to clean
- ▶ 15 % more microbeads to perfectly conform to the animal's shape
- ▶ The vacuum system allows complete immobilisation in the desired position
- ▶ Thermal insulation
- ▶ Safety valve can maintain vacuum up to 72 hours
- ▶ Warning: do not use with electric scalpel
- ▶ Supplied by unit with a pump

Small

4804964

Medium

4804965



GENIA COSYPAD®

Second generation meshed mat realized in a material used as anti-bedsore, Cosypad® offers comfort, hygiene and ease of cleaning. It constitutes an alternative of choice in duck-boards, underpads, and others litters for the wellbeing of animals. Thanks to its aerated structure, animals always stay clean and dry. Cosypad® satisfies all the conditions allying the comfort for our animals and the gain of time for professionals.

- ▶ Reduces the pressure exercised on muscles and articulations for an optimal comfort and the prevention of bedsores and urin scalds
- ▶ Hygienic, liquids are evacuated under the mat leaving the animal dry
- ▶ Hypoallergenic material, odorless, antibacterial
- ▶ Very resistant, rot-proof and insulating
- ▶ One non-slip surface
- ▶ Uses: in kennels, on tables, during post-operative wake up stage, for descaling (fast evacuation of liquids), as anti-fatigue mat (for the surgeon)
- ▶ Machine washable in front loaders only, to 40 °C



Size 1

- ▶ L 13" x W 19.7"

4809465

Size 3

- ▶ L 21.7" x W 28"

4809467

Size 2

- ▶ L 15.4" x W 28"

4809466

Size 4

- ▶ L 23.6" x W 45.3"

4809469

GENIA HANDY CANIS

Adjustable lifting straps and washable hand-held harness, for mild to serious temporary or permanent paresis or paralysis of the back legs. Helps the animal move about, especially climbing up stairs. Front leg Handy Canis also available.

Size XS

- ▶ Dimension A (in cm): 30–33
- ▶ Dimension B (in cm): 7.5–17
- ▶ Breed: Yorkshire, Chihuahua

4808526

Size S

- ▶ Dimension A (in cm): 33–36
- ▶ Dimension B (in cm): 15–23
- ▶ Breed: Dachshund

4808527

Size M

- ▶ Dimension A (in cm): 36–50
- ▶ Dimension B (in cm): 22–32
- ▶ Breed: Basset Hound

4808528

Size L


- ▶ Dimension A (in cm): 50–63
- ▶ Dimension B (in cm): 31–36
- ▶ Breed: Golden Retriever, German Shepherd

4808529

Size XL

- ▶ Dimension A (in cm): 63–81
- ▶ Dimension B (in cm): 37–53
- ▶ Breed: St. Bernard, Newfoundland

4808530



Sizing



Application



4808528

Handy Canis Front Legs

Size S

4808522

Size L

4808524

Size M

4808523

Size XL

4808525

EICKSKIN SKIN STAPLER

The EickSkin Disposable Skin Stapler is a simple and affordable stapler that provides consistent and reliable performance. Designed to fit comfortably in your hand, it enhances control of skin closure reducing surgery time as well as minimising infections and scarring. In addition the EickSkin Staples can be easily removed causing no pain to the patient.

- ▶ Pack of 10
- ▶ 35 x 7 mm staples
- ▶ Individually packaged – sterile

4801161



4801161

Skin Staple Remover

- ▶ Dedicated surgical staple remover

4803613

SURGICAL DRAPES

- ▶ Surgical grade, plain drape
- ▶ Durable cotton material
- ▶ Washable and autoclavable
- ▶ Colour: green

Dimensions (in cm): L 40 x W 60


4801007

Dimensions (in cm): L 70 x W 70

4801008

Dimensions (in cm): L 73 x W 90

4801009



4801007

SURGICAL DRAPE – FENESTRATED

- ▶ Surgical grade, fenestrated drape
- ▶ Durable cotton material
- ▶ Washable and autoclavable
- ▶ Colour: green

▶ Fenestration: 15 cm

▶ Dimensions (in cm): L 35 x W 60

4801010

▶ Fenestration: 20 cm


▶ Dimensions (in cm): L 73 x W 70

4801011

▶ Fenestration: 35 cm

▶ Dimensions (in cm): L 76 x W 90

4801012



4801010

EICKEMEYER® I.V. GIVING SET

- ▶ 20 drops/ml
- ▶ Length: 1.5 m
- ▶ Y-connector and Filter
- ▶ Includes Luer-Lock connector
- ▶ Individually wrapped sterile
- ▶ Box of 100
- ▶ DEHP & latex free


4802107



EICKEMEYER® I.V. GIVING SET EXTENSION

- ▶ Length: 1.5 m
- ▶ Individually wrapped sterile
- ▶ Pack of 50
- ▶ Sold individually

4804401



ANAESTHESIA MASKS FOR CATS & DOGS

- ▶ Robust, transparent plastic
- ▶ Size 3–6 incl. fixing trap

Size 1

- ▶ Ø opening: 22 mm
- ▶ Height: 24 mm
- ▶ Ø outside: 44 mm

4802594

Size 2

- ▶ Ø opening: 30 mm
- ▶ Height: 29 mm
- ▶ Ø outside: 57 mm

4802595



4802594



4802595

Size 3

- ▶ Ø opening: 33 mm
- ▶ Height: 35 mm
- ▶ Ø outside: 60 mm

4802596

Size 4

- ▶ Ø opening: 30 mm
- ▶ Height: 73 mm
- ▶ Ø outside: 87 mm

4802597



4802596



4802597

Size 5

- ▶ Ø opening: 42 mm
- ▶ Height: 97 mm
- ▶ Ø outside: 110 mm

4802598

Size 6

- ▶ Ø opening: 54 mm
- ▶ Height: 121 mm
- ▶ Ø outside: 130 mm

4802599



4802598



4802599

Complete Set (all 6 sizes)

4802593

REBREATHING BAGS

- ▶ Antistatic black rubber
- ▶ Includes: Hanging hook

0.5 Litre

4800227

2.0 Litre

4800231

0.75 Litre

4800228

3.0 Litre

4800232

1.0 Litre

4800229

4.0 Litre

4800233

1.5 Litre

4800230

5.0 Litre

4800234



4800227



4800231

PVC ENDOTRACHEAL TUBES

- ▶ Valve for Luer and Luer lock fittings
- ▶ Flattened, atraumatic tip, inflatable cuff
- ▶ Includes connector
- ▶ Cuff

- ▶ Ø inside: 2.5 mm
- ▶ Ø outside: 4 mm
- ▶ Length: 140 mm

4808343

- ▶ Ø inside: 3 mm
- ▶ Ø outside: 4.2 mm
- ▶ Length: 160 mm

4808344

- ▶ Ø inside: 3.5 mm
- ▶ Ø outside: 4.8 mm
- ▶ Length: 180 mm

4808345

- ▶ Ø inside: 4 mm
- ▶ Ø outside: 6 mm
- ▶ Length: 200 mm

4808116

- ▶ Ø inside: 4,5 mm
- ▶ Ø outside: 6.2 mm
- ▶ Length: 220 mm

4808117

- ▶ Ø inside: 5 mm
- ▶ Ø outside: 6.9 mm
- ▶ Length: 240 mm

4804941

- ▶ Ø inside: 5.5 mm
- ▶ Ø outside: 7.3 mm
- ▶ Length: 270 mm

4804942

- ▶ Ø inside: 6 mm
- ▶ Ø outside: 8 mm
- ▶ Length: 280 mm

4804943

- ▶ Ø inside: 6.5 mm
- ▶ Ø outside: 8,7 mm
- ▶ Length: 290 mm

4804944

- ▶ Ø inside: 7 mm
- ▶ Ø outside: 9.3 mm
- ▶ Length: 300 mm

4804945

- ▶ Ø inside: 7.5 mm
- ▶ Ø outside: 10 mm
- ▶ Length: 310 mm

4804946

- ▶ Ø inside: 8 mm
- ▶ Ø outside: 10.9 mm
- ▶ Length: 320 mm

4804947

- ▶ Ø inside: 8.5 mm
- ▶ Ø outside: 11.3 mm
- ▶ Length: 330 mm

4804948

- ▶ Ø inside: 9 mm
- ▶ Ø outside: 12 mm
- ▶ Length: 340 mm



4804949

- ▶ Ø inside: 9.5 mm
- ▶ Ø outside: 12.7 mm
- ▶ Length: 350 mm

4804950

- ▶ Ø inside: 10 mm
- ▶ Ø outside: 13.6 mm
- ▶ Length: 360 mm

4804951




4808343

INSTRUMENT MARKING TAPE SET

- ▶ For colour coding all types of instruments
- ▶ Autoclavable
- ▶ Dimensions: L 6 m x W 6.2 mm
- ▶ 8 different marking tapes

4801006



SANDBAGS

- ▶ Position aid for patients requiring x-rays or surgery
- ▶ Red, durable vinyl cover
- ▶ Malleable
- ▶ Easy to clean and disinfect
- ▶ *Please note: Not safe for washing machines

Small

- ▶ Weight: 1.2 kg
- ▶ Dimensions (in cm):
L 25 x W 15

4801663

Large

- ▶ Weight: 3.5 kg
- ▶ Dimensions (in cm):
L 30 x W 24

4801665

Medium

- ▶ Weight: 2 kg
- ▶ Dimensions (in cm):
L 27 x W 18

4801664



4801663

4801664

4801665


TOURNIQUETS

- ▶ Metal locking device
- ▶ Quick Release

4800131

- ▶ Plastic locking device
- ▶ Quick Release

4809888



4800131



4809888

PET NAIL CLIPPER

Pet Nail Clipper

- ▶ Large
- ▶ Spring action

4808333

Pet Nail Clipper Scissors

4808334



4808333



4808334

BITENOT COLLARS

Size 00

- ▶ Cats (< 4.5 kg)
- ▶ Dimensions (in mm):
L 360 x H 64

4801305

Size 0

- ▶ Cats (> 4.5 kg)
- ▶ Dimensions (in mm):
L 380 x H 76

4801306

Size 1

- ▶ Yorkshire Terrier (> 4.5 kg)
- ▶ Dimensions (in mm):
L 406 x H 89

4801307

Size 2

- ▶ Poodles (4.5 to 6.8 kg)
- ▶ Dimensions (in mm):
L 457 x H 102

4801308

Size 3

- ▶ Cocker Spaniels (8 to 15 kg)
- ▶ Dimensions (in mm):
L 534 x H 127

4801309

Size 3 1/2

- ▶ Spaniels (12 to 18 kg)
- ▶ Dimensions (in mm):
L 610 x H 127

4801313

Size 4

- ▶ Labrador (16 to 32 kg)
- ▶ Dimensions (in mm):
L 686 x H 152

4801310

Size 5

- ▶ Rottweilers (32 to 45 kg)
- ▶ Dimensions (in mm):
L 712 x H 178

4801311

Size 6

- ▶ St. Bernard (46 to 63 kg)
- ▶ Dimensions (in mm):
L 813 x H 203

4801312



Application



4801307

GENIA DIGIFLASH THERMOMETER

The Genia Digiflash Thermometer for rectal use, and is suitable for house pets, horses, cattle, sheep, pigs and other livestock. Supplies a reading within 5 seconds. This easy-to-use device is waterproof and has a flexible tip. Supplied with extra battery.

- ▶ Flexible tip
- ▶ Temperature measurement in 5 seconds
- ▶ Temperature range: 32 to 42 °C / 90 to 108 °F
- ▶ Accuracy: 0.1 °C / 0.18 °F
- ▶ Alarm indicating the end of measurement
- ▶ Low battery symbol display
- ▶ Auto power-off
- ▶ Waterproof, shock resistant
- ▶ Tip length: 75 mm
- ▶ Total length: 155 mm


4808531



GENIA PILL APPLICATOR

- ▶ PVC pill applicator with flexible tip
- ▶ Atraumatic soft tip
- ▶ Can be used with just one hand
- ▶ Add water if needed
- ▶ Supplied by 12 in a dispenser box

4809472



GENIA TABLET CRUSHING CUPS

- ▶ Plastic cup for cutting, crushing and storing tablets
- ▶ Ø 40 mm
- ▶ Height: 64 mm

4809473



GENIA K-STAR®

- Infusion tubing anti-pinch device
- Prevents pinching of infusion tubing
- Helps to secure the catheter
- Easy to use
- Quick to set up
- Dimensions (in cm): L 3 x W 3
- Pack of 50

4809060



Application



Application

GENIA NEOCOL COLLARS

- Transparent colourless polypropylene E-collars with quick fastening system and shock absorbing reflective neoprene outside edge
- Easy storage system
- 12 collars per pack

7.5 cm
4808040

20 cm
4808044

10 cm
4808041

25 cm
4808045


12.5 cm
4808042

30 cm
4808046

15 cm
4808043

Neocol Collar Set

- 2 each size
4808047



4808047

GENIA FUNCOL COLLARS

- ▶ These bright collars give owners the feeling of doing something special for their pets
- ▶ High quality collar with colourful paw prints
- ▶ Easy to put on pets with quick fastening system
- ▶ Easy storage system
- ▶ 12 collars per pack

7.5 cm, pink

4808053

10 cm, purple

4808054

12.5 cm, orange

4808055

15 cm, green

4808056

20 cm, yellow

4808057

25 cm, red

4808058

30 cm, blue

4808059

Funcol Collar Set

- ▶ 2 each size

4808060



4808060

GENIA GRENCOL COLLARS

- ▶ Transparent colourless polypropylene E-collars with quick fastening system
- ▶ Wavy outer edge avoids sharp frontal impact
- ▶ Easy storage system
- ▶ Made from 80 % recycled material!
- ▶ 12 collars per pack

7.5 cm

4808512

10 cm

4808513

12.5 cm

4808514

15 cm

4808515

20 cm

4808516

25 cm

4808517

30 cm

4808518

35 cm

4808520

40 cm

4808521

Grencol Collar Set

- ▶ 2 each size
(except 35/40 cm collar)

4808519



4808519

I.V. DRIP STAND AND ACCESSORIES

I.V. Drip Stand

- ▶ Chrome finish
- ▶ 5 Star base
- ▶ 4 Prong hooks
- ▶ 5 Swivelling castors
- ▶ 2 Lockable castors
- ▶ Adjustable height: 120 – 210 cm

4800987

Storage shelf

- ▶ For I.V. Drip Stand (4800987)
- ▶ For pipe diameters up to 19 mm
- ▶ Dimensions (in cm): L 20 x W 15

4808822



4800987



4808822

CANINE SPAY KIT

Ovariohysterectomy is one of the most common procedures performed in veterinary clinics. This procedure helps prevent mammary tumours, pyometra and unwanted pregnancy. This kit provides all the necessary equipment for veterinary surgeons to perform a bitch spay.


*All of these instruments may be purchased separately.

*Large assortment of surgical instruments available.

Complete set, consisting of:

- ▶ Surgical Scalpel Handle No. 3, blade sizes 10 to 15, L 125 mm
- ▶ OLSEN-HEGAR Needle Holder (TC), EICKEMEYER® Lifetime Warranty, Tungsten Carbide, L 140 mm
- ▶ NELSON-METZENBAUM Scissors, curved, L 145 mm
- ▶ MAYO Scissors, curved, L 170 mm
- ▶ Dressing Forceps, standard tip, L 160 mm
- ▶ ADSON-BROWN Tissue Forceps, delicate, 9 x 9 teeth, L 120 mm
- ▶ 2x DEBAKEY-MOSQUITO Artery Forceps, curved, L 120 mm
- ▶ 2x KELLY Artery Forceps – Basic, curved, L 140 mm
- ▶ ROCHESTER-CARMALT Artery Forceps – Basic, straight, L 160 mm
- ▶ ROCHESTER-CARMALT Artery Forceps – Basic, curved, L 160 mm
- ▶ 2x ALLIS Tissue Forceps, 5 x 6 teeth, L 150 mm
- ▶ 4x BACKHAUS Towel Clip, L 110 mm
- ▶ Spay Hook, Ø button tip: 3.0 mm, L 160 mm
- ▶ Premium Instrument Box, made of aluminium, perforated bottom and lid, features side closures, not autoclavable, dimensions (cm): L 30 x W 14 x H 7

4800193



AESCULAP ISIS CAT CLIPPER

The AESCULAP Isis cordless clipper is ideal for preparation for small areas (forelimb, head and neck) when taking bloods, biopsy or for hard-to-reach areas.

- ▶ Quiet and low vibration
- ▶ Ergonomic design
- ▶ Cutting width: 24 mm
- ▶ Cutting length: 0.5 mm
- ▶ Li-ion battery: runtime up to 120 min with approx. 90 min charging time
- ▶ Charge indicator via charging station
- ▶ Weight: 120 g (incl. battery & blade)
- ▶ Includes: clipper, blade (4802265), (carbon blade – cannot be sharpened), charging station, attachment comb (adjustable to 4 different cutting lengths, 3–6 mm), brush, AESCULAP oil, instructions

4802162



AESCULAP VEGA DOG CLIPPER

Superior cordless clipping. Reduced noise and vibration, whilst delivering superior cordless clipping. Ideal for abdominal, forelimbs and neck regions.

- ▶ Blade cutting width: 40 mm
- ▶ Adjustable blade cutting length: 0.5/1.5/2.5 mm
- ▶ Battery life: Up to 50 min/2 hour recharge time
- ▶ Charging indicator
- ▶ Weight: 250 g
- ▶ Includes: Vega clipper blade (4802266), charging station

4802164



ZWEEZI SADDLE STOOL

- ▶ Padded synthetic leather seat
- ▶ Height adjustable: 55–85 cm
- ▶ Gas lift mechanism
- ▶ Base diameter: 54 cm
- ▶ 5 star base with swivable castors
- ▶ Colour: evergreen

4800594

